APTIC[®] SURFACE PLUS

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APTIC[®] SURFACE PLUS

High-performance concentrated surface disinfectant

For Operating Rooms and Intensive Care Units

APTIC[®] Surface Plus is a next-generation concentrated disinfectant based on Quaternary Ammonium Compounds (QACs). Specifically designed for surface disinfection in high-sensitivity environments such as Operating Rooms (ORs), Intensive Care Units (ICUs), Burn Units, and Isolation Wards. Its high concentration allows economic dilution according to the required level of contamination, maintaining excellent efficacy against bacteria, viruses, and fungi

Technical Features

- ✓ Proven efficacy against: MRSA, VRE, Pseudomonas aeruginosa, Escherichia coli, Staphylococcus aureus, Candida albicans, Influenza virus, SARS-CoV-2 (COVID-19).
- ✓ High concentration - Economical dilution ratios.
- ✓ Short contact time (5-10 minutes depending on dilution).
- ✓ Safe on surfaces - Non-corrosive to metals, plastics, glass, and paints.
- ✓ Chlorine-free, alcohol-free, glutaraldehyde-free - No unpleasant odor or respiratory irritation.
- ✓ Effective in the presence of organic load (blood, protein, secretions).
- ✓ Leaves antimicrobial residue that prevents recontamination for several hours.

Key Applications

Area	المجال
Operating Rooms	Operating tables, anesthesia equipment, instrument surfaces, floors, walls
Intensive Care Units (ICU)	Patient beds, ventilators, monitoring screens, service tables
Burn & Isolation Units	Patient surrounding surfaces, room floors, door handles
Medical Laboratories	Work benchtops, incubators, centrifuges
Outpatient Clinics	Examination surfaces, dental chairs, waiting area floors

Directions for Use (Dilutions):

Application	Dilution Ratio	Contact Time	Notes
Critical surface disinfection (OR)	(ml per liter of water 10) 1:100	10 minutes	Wipe and leave to dry
Routine surface disinfection	(ml per liter of water 5) 1:200	5 minutes	For non-critical surfaces
Heavily soiled surfaces	(ml per liter of water 20) 1:50	15 minutes	For blood or heavy secretions
Floors and walls	(ml per liter of water 5) 1:200	5 minutes	Mop and leave to air dry

Safety & Storage

- ✓ Store in original closed container at 5-30°C
- ✓ Avoid contact with eyes - if contacted, rinse generously with water for 15 minutes.
- ✓ Use appropriate gloves in critical environments.
- ✓ Do not mix with anionic soaps or strong alkaline detergents.
- ✓ Diluted solution is stable for 7 days.